



Rx Only

Caution: Federal law (USA) restricts this product to sale by or on the order of a physician (or properly licensed healthcare professional).

OCM™

INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION

OCM™ is a wound care drug-device comprised of fish peptides infused with cod liver oil, which acts as an anhydrous skin protectant. When applied to a wound surface, the matrix is naturally incorporated into the wound over time. OCM™ is designed for intimate contact with both regular and irregular wound beds, to provide a conducive environment for the patient's natural wound healing process.

INDICATIONS FOR USE

OCM™ is indicated for the management of wounds including: partial and full-thickness wounds, pressure ulcers, venous ulcers, diabetic ulcers, chronic vascular ulcers, tunneled/undermined wounds, surgical wounds (donor sites/grafts, post-Moh's surgery, post-laser surgery, podiatric, wound dehiscence), trauma wounds (abrasions, lacerations, superficial partial thickness burns, skin tears) and draining wounds.

CONTRAINDICATIONS

- This product should not be used in patients with known sensitivity to fish products or cod liver oil.
- The product is not indicated for use in third degree burns.

PRECAUTIONS

- Each unit of OCM™ is for single use. Discard all opened and unused portions after each treatment session.
- Contents are sterile if the inner package is unopened and undamaged. Do not use if the inner package is broken. Do not re-sterilize.
- Discard product if mishandling has caused possible damage or contamination.
- Excessive exudate, bleeding, acute swelling, and infection should be controlled before OCM™ is applied. Debridement or excision must be done thoroughly to remove any remaining necrotic tissue that may cause infection.
- The following complications are possible with the use of wound management products: infection, chronic inflammation (initial application of wound dressings may be associated with transient, mild, localized inflammation), allergic reaction, excessive redness, pain, or swelling. If any of these conditions occur, the product should be removed.

DIRECTIONS FOR USE

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

- Up to 3.2g OCM™ may be applied daily.
- Up to 6.4g OCM™ may be applied weekly.

Initial application

- Always handle OCM™ using aseptic techniques. Only inner package content is sterile.
- Prior to application of OCM™, prepare the wound bed using standard methods to ensure wound is free of debris and necrotic tissue. If necessary, surgically debride the wound to ensure the wound edges contain viable tissue.
- Snap open the tip of the OCM™ package.
- Slowly squeeze to apply either directly to the wound bed in stripes, or onto a suitable applicator (e.g. a tongue depressor). Distribute product evenly across the wound bed so that it is uniformly applied and fills the wound bed. Cover the periwound to ensure complete coverage.
- After application, cover with an appropriate dressing to maintain OCM™ adherence and protect the wound area. The optimum dressing is determined by wound location, size, depth, and volume of exudate.

Follow-up application

- Duration of treatment and reapplication frequency is determined by the physician and will depend on the condition of the patient as well as the level of wound exudate.
- Reapply when OCM™ has been naturally incorporated into the wound bed or as directed by a wound care professional.
- Remove the dressing from the wound bed gently.
- Cleanse the wound area and prepare the wound bed prior to application.
- Distribute product evenly across the wound bed to fill the wound and cover the periwound to ensure complete coverage.
- Cover with an appropriate dressing to maintain OCM™ adherence and protect the wound area.
- Change the dressing based on physician assessment of patient's individual needs and as consistent with manufacturer's instructions to maintain a moist, clean wound area.

HOW SUPPLIED

| Product Codes | Size | Quantity |
|---------------|-------|----------------|
| CM-720.001 | 1.6 g | 1 device/pouch |

STORAGE

Store at room temperature [77°F/25°C]. Keep away from sunlight. See product pouch for expiration date.

SYMBOLS USED IN LABELING

| ISO 15223-1: 2016 Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements | | |
|---|-------|----------------------------------|
| | 5.1.6 | Catalogue Number |
| | 5.1.5 | Batch Code |
| | 5.4.2 | Do not re-use |
| | 5.2.6 | Do not re-sterilize |
| | 5.2.8 | Do not use if package is damaged |
| | 5.1.4 | Use-by date |
| | 5.2.4 | Sterilized using irradiation |
| | 5.4.3 | Consult instructions for use |
| | 5.3.2 | Keep away from sunlight |
| | 5.4.4 | Caution |
| | 5.1.1 | Manufacturer |
| | 5.3.6 | Upper limit of temperature |

For product ordering information, further information, or guidance, please contact Omeza for assistance at: 1-888-88OMEZA [1-888-886-6392]

PRODUCT INFORMATION DISCLOSURE

OMEZA LLC HAS EXERCISED REASONABLE CARE IN THE SELECTION OF MATERIALS AND THE MANUFACTURE OF THESE PRODUCTS. OMEZA LLC EXCLUDES ALL WARRANTIES, WHETHER EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. OMEZA LLC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM USE OF THIS PRODUCT. OMEZA LLC NEITHER ASSUMES NOR AUTHORIZES ANY PERSON TO ASSUME FOR IT ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THESE PRODUCTS.

Omeza LLC
1610 Northgate Blvd
Sarasota, FL 34234

www.omeza.com

Omeza® is a registered trademark of Omeza LLC in the United States and/or other countries. © 2024 Omeza LLC.